

VELP Scientifica MSL Digital Series



MSL 8 Digital



MSL 25 Digital



MSL 50 Digital

Stirring speed rpm	40 to 1500	40 to 1500	20 to 1100
Stirring volume L (H ₂ O)	up to 8	up to 25	up to 50
Timer	•	•	•
Reverse Operation	•	•	•
Intermittent Mode	•	•	•
Dimensions (WxHxD) mm	168x60x183	253x71x269	354x77x369
Weight Kg	0.95	1.8	4
Power Supply	100÷240 V	100÷240 V	100÷240 V
Power input	5 W	9 W	10 W



www.velp.com