

## VELP Scientifica Solutions for BOD Determination



**BOD Sensor Set**



**BOD Sensor System 6**



**BOD Sensor System 10**



**BOD EVO Sensor System 6 \***

Real Time Data Reading (ppm BOD)		On Display	On Display	On Display	On Display
PC Data Reading (ppm BOD)		-	-	-	with Wireless DataBox™ and BODSoft™
Data Storage	for a 5-day test (BOD <sub>5</sub> )	Every 24 hours	Every 24 hours	Every 24 hours	Every 30 min/1/2/4/6/8/12/24 hours
	for a longer test	Last 5 values	Last 5 values	Last 5 values	Every 2/4/6/8/12/24 hours
Scales (ppm BOD)		90, 250, 600, 999 ppm	90, 250, 600, 999 ppm	90, 250, 600, 999 ppm	90, 250, 600, 999 ppm
Supplied with	Sensor	1	6	10	6 EVO
	Dark Glass Bottle	1	6	10	6
	Alkali Holder	1	6	10	6
	Stirring Bar	1	6	10	6
	Stirring Station	-	6-position	10-position	6-position
Dimensions (WxHxD)	mm	75x230x75	270x300x185	432x300x165	270x300x185
	(in)	(3.0x9.3x3.0)	(10.6x11.8x7.3)	(17.0x11.8x6.0)	(10.6x11.8x7.3)
Weight	Kg	0.4	3.4	6.5	3.4
	(lb)	(0.9)	(7.5)	(14.5)	(7.5)
Power Supply		-	115 or 230 V	115 or 230 V	115 or 230 V
Power		2 lithium batteries	2 W	2 W	2 W

\* Wireless DataBox™ and BODSoft™ are included for those codes that end with 46.



[www.velp.com](http://www.velp.com)